


ENDOCOR
PIONEERING INNOVATION FOR LIFE™

HYPERFLOW™

Compact & Ergonomic Inflation Device



30 ATM PRESSURE GAUGE
20 CC SYRINGE
90° ROTARY PRESSURE GAUGE

HYPERFLOW™

Compact & Ergonomic Inflation Device

PRESSURE GAUGE

- 30 atm
+ reliable in any procedure
- 90° rotating manometer
+ "custom" handling
- coded figures from nominal to burst pressures
+ easy reading
- luminescent dial
+ great visibility

BODY & BUTTON

- unique mechanism for rapid pressure increase
+ efficiency, fast procedures
- 20 cc graduated syringe
+ compact system
- grip design
+ better handling
- clear material
+ easy air bubble detection and removal

HANDLE

- ergonomic design
+ optimal handling
- progressive resistance while inflating
+ pressure feedback

Manufacturer:
ENDOCOR GmbH
GERMANY
info@endocor.com
www.endocor.com

PRODUCT ORDERING INFORMATION

HYPERFLOW™	Accessories	Article-#	Packaging	Shelf life
Inflation Device	3-way stopcock	HF3020	12-unit box	5 years

CE 0481

© ENDOCOR

All rights reserved. Technical specification and data are subject to modification, revision and improvement without notice. RJIS20012016