

**ENDOCOR**  
PIONEERING INNOVATION FOR LIFE™

# ENDOGLIDE DF **Plus**

PTA Dialysis Fistula Balloon Catheter Plus



**FOR SPECIAL CASES**  
OF DIALYSIS FISTULA WITH  
HIGH PRESSURE NEEDS

# ENDOGLIDE DF Plus™

An exceptional & highly specialized balloon for the most difficult cases where additional pressure is a must.

INCREASED PUSHABILITY FOR A SIGNIFICANT AMOUNT OF APPLIED FORCE TRANSITION FROM THE HUB TO THE TIP



## HIGHLIGHTS

- Super low entry profile
- OTW System with 0.018" or 0.035" guidewire compatibility
- Very high pressure with rated burst pressure of 26 bar
- Wide working pressure range between 10 and 35 bar (ABP)
- High visibility with 2 Gold marker bands
- Smooth and safe handling through precision laser-bond at critical flex points
- Maximum flexibility for easy navigation reaches the most distal lesions
- High Rated Burst Pressure optimizes working range and increases the safety level for patients

# ENDOGLIDE DF Plus™


PTA Dialysis Fistula Balloon Catheter Plus

## PRODUCT ORDERING INFORMATION

PRODUCT RANGE FOR ENDOGLIDE DF Plus™				
Usable catheter length 45 cm				
Balloon diameter (mm)	Balloon length (mm)			
	20	30	40	60
4	45EH420	45EH430	45EH440	45EH460
5	45EH520	45EH530	45EH540	45EH560
6	45EH620	45EH630	45EH640	45EH660
7	45EH720	45EH730	45EH740	45EH760
8	45EH820	45EH830	45EH840	45EH860
9	45EH920	45EH930	45EH940	45EH960
10	45EH1020	45EH1030	45EH1040	45EH1060

## TECHNICAL SPECIFICATIONS

Mode of exchange	OTW	Usable catheter length	45 cm
Balloon characteristic	semi-compliant	Nominal pressure	10 bar
Balloon material	modified polyamide	Rated burst pressure	26 bar
Guidewire compatibility (max)	0.035"	Crossing profile	ø 4.0 - 5.0 < 2.00 mm ø 6.0 - 7.0 < 2.10 mm ø 8.0 - 10.0 < 2.20 mm
Sheath compatibility (min)	5F (ø 4.0 - 5.0 mm)	Average burst pressure	35 bar
	6F (ø 6.0 - 7.0 mm)	Total shelf life	4 years
	7F (ø 8.0 - 10.0 mm)		



# ENDOGLIDE DF **Plus**™

FOR SPECIAL CASES  
OF DIALYSIS FISTULA WITH  
HIGH PRESSURE NEEDS



Manufacturer:  
**ENDOCOR** GmbH  
GERMANY  
[info@endocor.com](mailto:info@endocor.com)  
[www.endocor.com](http://www.endocor.com)

Distributed by:



**CE** 0481

© ENDOCOR

All rights reserved. Technical specification and data are subject to modification, revision and improvement without notice. RJIS20012016