

ENDOCOR
PIONEERING INNOVATION FOR LIFE™

ENDOGLIDE DF

PTA Dialysis Fistula Balloon Catheter




ENHANCED SAFETY &
LONG TERM OUTCOMES
FOR DIALYSIS PATIENTS

ENDOGLIDE DF™

Designed for patients who depend on frequent dialysis treatment, the ENDOGLIDE DF facilitates vascular access and ensures safer and shorter treatment of stenosis and fistula thrombosis.

SPECIFICALLY DESIGNED TO SUCCESSFULLY DILATE EVEN THE MOST RESILENT STENOSIS THAT ARE NORMALLY RESISTANT TO EXPANSION



HIGHLIGHTS

- Semi-compliant Over-the-Wire (OTW) balloon catheter
- High pressure 18 bar
- High rated burst pressure optimizes working range and increases the safety level for patients
- Fast inflation and deflation time maximizes efficiency during treatment and reduces the procedure time compared to other industry participants
- Reduced radiation exposure for clinician, team and patient due to decreased time for balloon operation
- Lowest entry profile in the industry with the addition of a beveled tip allows access of the most difficult lesions
- Durable balloon material and shaft with high pushability
- Maximum flexibility for easy navigation reaches the most distal lesions

ENDOGLIDE DF™

PTA Dialysis Fistula Balloon Catheter 0.018" HP

PRODUCT ORDERING INFORMATION - UCL 45 CM

4 mm Diameter	
Length (mm)	Catalogue Number
20	04EE040020BTO
40	04EE040040BTO

5 mm Diameter	
Length (mm)	Catalogue Number
20	04EE050020BTO
40	04EE050040BTO

6 mm Diameter	
Length (mm)	Catalogue Number
20	04EE060020BTO
40	04EE060040BTO

7 mm Diameter	
Length (mm)	Catalogue Number
20	04EE070020BTO
40	04EE070040BTO

PRODUCT ORDERING INFORMATION - UCL 80 CM

4 mm Diameter	
Length (mm)	Catalogue Number
20	08EE040020BTO
40	08EE040040BTO


5 mm Diameter	
Length (mm)	Catalogue Number
20	08EE050020BTO
40	08EE050040BTO

6 mm Diameter	
Length (mm)	Catalogue Number
20	08EE060020BTO
40	08EE060040BTO

7 mm Diameter	
Length (mm)	Catalogue Number
20	08EE070020BTO
40	08EE070040BTO

TECHNICAL SPECIFICATIONS


Guidewire	0.018"
Mode of exchange	over-the-wire (OTW)
Usabel catheter length	45 cm, 80 cm
Nominal pressure	9 bar
Rated burst pressure	18 bar
Coating	hydrophilic coating
Balloon material	nylon
Balloon characteristic	semi compliant
Platinum X-ray balloon marker	2 marker - at proximal and distan end of balloon
Low entry profile	0.016" (0.40 mm)



FLEXIBLE EVEN AT VERY HIGH PRESSURE

ENDOGLIDE DF

ENHANCED SAFETY &
LONG TERM OUTCOMES
FOR DIALYSIS PATIENTS



Manufacturer:
ENDOCOR GmbH
GERMANY
info@endocor.com
www.endocor.com

Distributed by:



CE 0481

© ENDOCOR

All rights reserved. Technical specification and data are subject to modification, revision and improvement without notice. RJIS20012016