


ENDOCOR
PIONEERING INNOVATION FOR LIFE™

ENDOGLIDE HP CTO


PTA High Pressure CTO Balloon Catheter



THE CHOICE FOR
EFFECTIVE TREATMENT
OF HIGHLY RESILIENT STENOSIS

ENDOGLIDE HP CTO™

DESIGNED TO PERFORM RELIABLY IN BOTH EVERY DAY & WORK HORSE CASES AS WELL AS THE MOST CHALLENGING OF SITUATIONS.



HIGHLIGHTS

ENDOGLIDE HP CTO™ IS A PTA HIGH PRESSURE BALLOON CATHETER SPECIFICALLY DESIGNED FOR DILATATION IN COMPLEX CASES. THIS DEVICE EFFECTIVELY HELPS TO TREAT SEVERE STENOSIS THAT ARE FREQUENTLY RESISTANT TO EXPANSION.

- OTW System with 0.035" guidewire compatibility
- High pressure with rated burst pressure of 26 bar
- Wide working pressure range between 10 and 35 bar
- High visibility with 2 Gold marker bands
- Smooth and safe handling through precision laser-bond at critical flex points

ENDOGLIDE HP CTO™

PTA High Pressure CTO Balloon Catheter


PRODUCT RANGE FOR ENDOGLIDE HP CTO™				
Usable catheter length 45 cm				
Balloon diameter (mm)	Balloon length (mm)			
	20	30	40	60
4	45EH420	45EH430	45EH440	45EH460
5	45EH520	45EH530	45EH540	45EH560
6	45EH620	45EH630	45EH640	45EH660
7	45EH720	45EH730	45EH740	45EH760
8	45EH820	45EH830	45EH840	45EH860
9	45EH920	45EH930	45EH940	45EH960
10	45EH1020	45EH1030	45EH1040	45EH1060

PRODUCT RANGE FOR ENDOGLIDE HP CTO™				
Usable catheter length 80 cm				
Balloon diameter (mm)	Balloon length (mm)			
	20	30	40	60
4	80EH420	80EH430	80EH440	80EH460
5	80EH520	80EH530	80EH540	80EH560
6	80EH620	80EH630	80EH640	80EH660
7	80EH720	80EH730	80EH740	80EH760
8	80EH820	80EH830	80EH840	80EH860
9	80EH920	80EH930	80EH940	80EH960
10	80EH1020	80EH1030	80EH1040	80EH1060

TECHNICAL SPECIFICATIONS


Mode of exchange	OTW
Balloon characteristic	semi-compliant
Balloon material	modified polyamide
Guidewire compatibility (max)	0.035"
Sheath compatibility (min)	6F (Ø 4.0 - 7.0 mm) 7F (Ø 8.0 - 9.0 mm) 8F (Ø 10.0 mm)

Usable catheter length	45 cm / 80 cm
Nominal pressure	10 bar
Rated burst pressure	26 bar
Crossing profile	ø 4.0 - 5.0 < 2.00 mm ø 6.0 - 7.0 < 2.10 mm ø 8.0 - 10.0 < 2.20 mm
Average burst pressure	35 bar
Total shelf life	4 years




ENDOGLIDE HP CTO

TRAVERSES HARD TO PASS
CHRONIC TOTAL OCCLUSIONS



Manufacturer:
ENDOCOR GmbH
GERMANY
info@endocor.com
www.endocor.com

Distributed by:



© ENDOCOR
All rights reserved. Technical specification and data
are subject to modification, revision and improvement
without notice. RJIS28062016