

ENDOCOR
PIONEERING INNOVATION FOR LIFE™

GI-GHOST PV™

BALLOON-EXPANDABLE PERIPHERAL STENT SYSTEM


SUPERIOR IN ITS CLASS
EXCELLENT FLEXIBILITY, PRECISE
PLACEMENT AND PREMIUM VISIBILITY

GHOST PV™

BALLOON-EXPANDABLE PERIPHERAL STENT SYSTEM

THIS BALLOON-EXPANDABLE PERIPHERAL STENT SYSTEM IS ONE OF ENDOCOR'S NEW GENERATION PRODUCTS.

IT IS MADE OF COBALT-CHROMIUM ALLOY AND PROVIDES A MEANS OF SAFELY ADVANCING THE STENT THROUGH THE PERIPHERAL VESSELS (E.G. EXTERNAL ILIAC, INTERNAL ILIAC, COMMON ILIAC, SUPERFICIAL FEMORAL, DEEP FEMORAL PROFUNDA, RENAL) TO THE DESIRED LOCATION.





HIGHLIGHTS



- LOW CROSSING PROFILE FOR ALL DIAMETERS
- COMBINATION LOW SHAFT AND TIP ENTRY PROFILE
- SOFT, TAPERED BALLOON TIP FOR LESS TRAUMA
- FLEXIBLE STENT DESIGN
- HIGH RADIAL STRENGTH
- EASY DEPLOYMENT AND PRECISE PLACEMENT
- EVEN AND ACCURATE STENT EXPANSION
- GOOD ACCESSIBILITY AND HIGH PUSHABILITY ACROSS THE PERIPHERAL VESSELS

BEATS ITS COMPETITORS WITH HIGHEST RATED BURST PRESSURE

BETTER THAN CORDIS & BOSTON SCIENTIFIC



BETTER THAN COOK & COVIDIEN



COBALT-CHROMIUM IS APPROXIMATELY 70-80% STRONGER THAN STAINLESS STEEL

GHOST PV	SF	SF
MODE OF EXCHANGE	Over-The-Wire (OTW)	
STENT PLATFORM	Cobalt-Chromium	
STENT LENGTHS	17 / 26 / 36 / 57 mm	
STENT DIAMETER	ø 5 / 6 / 7 / 8 / 9 / 10 mm	
STENT MATERIAL THICKNESS	145 µm	
DISTAL SECTION	Two coaxial lumen	
RECOMMENDED GUIDEWIRE	0.018" (0.46 mm)	0.035" (0.89 mm)
RECOMMENDED INTRODUCER SHEATH	5F	6F
BALLOON CHARACTERISTIC	Semi compliant	
BALLOON MATERIAL	Pebax 7233 SA01 / Nylon 12, transparent	
BALLOON FOLDING	3 folds	
NOMINAL PRESSURE	9 bar	
RATED BURST PRESSURE (RBP)	ø 5.0 to 7.0 mm ø 8.00 mm ø 9.0 and 10.0 mm	18 bar 16 bar 13 bar
USABLE CATHETER LENGTHS	80 / 135 cm	80 / 120 / 150 cm
X-RAY MARKER	2 Platinum- Iridium markers at each balloon end, proximal and distal	
ENTRY PROFILE	0.50 ± 0.1 mm (0.020")	0.95 ± 0.025 mm (0.037")
CROSSING PROFILE	max 1.87 mm (0.074")	max 2.00 mm (0.079")
HYDROPHILIC COATING	400 mm ± 15 mm	
SHELF LIFE	2 years	

BENEFIT:
GHOST PV IS MADE OF
COBALT-CHROMIUM

COMPETITOR
PRODUCTS MADE
OF STAINLESS STEEL

COMPANY	PRODUCT
BIOTRONIK	DYNAMIC
BARD	VALEO
COVIDIEN	PARAMOUNT MINI GPS™
CORDIS	PALMAZ GENESIS®
BOSTON SCIENTIFIC	EXPRESS SD®
COOK	FORMULA 418®





PRODUCT RANGE FOR **GIHOST PV™**

5F - UCL 80 cm

Stent Diameter in mm	Stent Length in mm			
	17	26	36	57
5.00	5GP17-80B	5GP26-80B	5GP36-80B	5GP57-80B
6.00	6GP17-80B	6GP26-80B	6GP36-80B	6GP57-80B
7.00	7GP17-80B	7GP26-80B	7GP36-80B	7GP57-80B
8.00	8GP17-80B	8GP26-80B	8GP36-80B	8GP57-80B
9.00	9GP17-80B	9GP26-80B	9GP36-80B	9GP57-80B
10.00	10GP17-80B	10GP26-80B	10GP36-80B	10GP57-80B

5F - UCL 135 cm

Stent Diameter in mm	Stent Length in mm			
	17	26	36	57
5.00	5GP17-13B	5GP26-13B	5GP36-13B	5GP57-13B
6.00	6GP17-13B	6GP26-13B	6GP36-13B	6GP57-13B
7.00	7GP17-13B	7GP26-13B	7GP36-13B	7GP57-13B
8.00	8GP17-13B	8GP26-13B	8GP36-13B	8GP57-13B
9.00	9GP17-13B	9GP26-13B	9GP36-13B	9GP57-13B
10.00	10GP17-13B	10GP26-13B	10GP36-13B	10GP57-13B



PRODUCT RANGE FOR **GIHOST PV™**

6F - UCL 80 cm

Stent Diameter in mm	Stent Length in mm			
	17	26	36	57
5.00	5GP17-80C	5GP26-80C	5GP36-80C	5GP57-80C
6.00	6GP17-80C	6GP26-80C	6GP36-80C	6GP57-80C
7.00	7GP17-80C	7GP26-80C	7GP36-80C	7GP57-80C
8.00	8GP17-80C	8GP26-80C	8GP36-80C	8GP57-80C
9.00	9GP17-80C	9GP26-80C	9GP36-80C	9GP57-80C
10.00	10GP17-80C	10GP26-80C	10GP36-80C	10GP57-80C

6F - UCL 120 cm

Stent Diameter in mm	Stent Length in mm			
	17	26	36	57
5.00	5GP17-12C	5GP26-12C	5GP36-12C	5GP57-12C
6.00	6GP17-12C	6GP26-12C	6GP36-12C	6GP57-12C
7.00	7GP17-12C	7GP26-12C	7GP36-12C	7GP57-12C
8.00	8GP17-12C	8GP26-12C	8GP36-12C	8GP57-12C
9.00	9GP17-12C	9GP26-12C	9GP36-12C	9GP57-12C
10.00	10GP17-12C	10GP26-12C	10GP36-12C	10GP57-12C


6F - UCL 150 cm

Stent Diameter in mm	Stent Length in mm			
	17	26	36	57
5.00	5GP17-15C	5GP26-15C	5GP36-15C	5GP57-15C
6.00	6GP17-15C	6GP26-15C	6GP36-15C	6GP57-15C
7.00	7GP17-15C	7GP26-15C	7GP36-15C	7GP57-15C
8.00	8GP17-15C	8GP26-15C	8GP36-15C	8GP57-15C
9.00	9GP17-15C	9GP26-15C	9GP36-15C	9GP57-15C
10.00	10GP17-15C	10GP26-15C	10GP36-15C	10GP57-15C

Distributed by:

Manufacturer:

ENDOCOR GmbH - GERMANY - info@endocor.com - www.endocor.com



© **ENDOCOR**

All rights reserved. Technical specification and data are subject to modification, revision and improvement without notice. RJIS02042018